



KeyGuard®

KeyGuard precision dental resin produces custom-fit mouthguards for unbeatable comfort and protection!



BIOCOMPATIBLE



STRONG



DURABLE



COMFORTABLE



WHITE

#4220200 (1kg bottle)



BLACK

#4220201 (1kg bottle)



PURPLE

#4220202 (1kg bottle)



GREEN

#4220203 (1kg bottle)



RED

#4220204 (1kg bottle)

FEATURES

- **BIOCOMPATIBLE**
- **PRECISE FIT FOR ENHANCED COMFORT AND PERFORMANCE**
- **PROTECTION ONLY WHERE NEEDED WITH COMPLETE DESIGN FREEDOM**
- **OPTIMIZED COST EFFICIENCY**

PROPERTY	METHOD	KEYGUARD®
Tear Strength	ASTM D624-00	30 kN/m
Hardness	ASTM D2240-00, 25C @ 0 sec	75-90 A
Impact Strength	ISO 179-1	27.5 kJ/m ²
Impact Resistance	ISO 179-1	220 J/m
Water Sorption	ASTM D570 @ 3hr	5.7%
Tensile Cyclic Hysteresis	Energy returned @ 25C Energy lost @ 25C	27% 73%
Elongation at break (strain)	ASTM D638	130%
Liquid Viscosity	Internal	900cP
Biocompatibility	ISO 10993-1	Pass

keystone
industries



keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

REV 0 05/2025