



KeyMask[®]


Flexible yet firm gingival mask that simulates the look and feel of gum tissue. Incorporate KeyMask into your digital restoration planning



ACCURATE



**LOW
SHRINKAGE**



FLEXIBLE



**SIMULATES
GUM TISSUE**

#4220002 (0.5kg bottle)

PROPERTY	METHOD	KEYMASK
Tensile Strength	ASTM D638	2 MPa
Tensile Young's Modulus	ASTM D638	4.6 MPa
Elongation at Break (%)	ASTM D638	41%



FLEXIBLE YET FIRM

This gingival mask simulates the look and feel of gum tissue



EASILY TRIMMED

Creates a proper emergence profile for pontic shape



keystone
industries

keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027

EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131

EU: +49 77 31 91 21 01

REV 0 05/2025