





KeyModel Ultra®


KeyModel Ultra is a material designed for 3D printing of dental and orthodontic models.



**ULTRA-FAST
TECHNOLOGY**



**FLAWLESS
SURFACE DETAIL**



**THERMOFORMING
RELEASE AGENT**



TPO FREE



 Ivory	1kg.....4220009	 Light Grey	1kg.....4220010
	4kg.....4240009		4kg.....4240010

PROPERTY	METHOD	KEYMODEL ULTRA
Maximum Tensile Strength	ASTM D638	50 MPa
Tensile Modulus	ASTM D638	1700 MPa
Elongation at Break	ASTM D638	5%
Maximum Flexural Strength	ASTM D790	70 MPa
Flexural Modulus	ASTM D790	1940 MPa
Viscosity @ 25C	ASTM D7867	500-600 cP

keystone
industries



keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027

EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131

EU: +49 77 31 91 21 01

REV 0 05/2025