



KeyModel[®]


KeyModel produces highly precise and accurate dental models. Its low shrinkage formulation makes it an ideal resin for crown, bridge, and implant models, orthodontic clear aligner models, and diagnostic modeling.



ACCURATE



**LOW
SHRINKAGE**



EASY VISIBILITY

#4200006 (0.5kg bottle)

#4220006 (1kg bottle)

PROPERTY	METHOD	KEYMODEL
Flexural Strength	ASTM D790	79.6 MPa
Flexural Modulus	ASTM D790	2293 MPa
Tensile Strength	ASTM D638	51
Tensile Young's Modulus	ASTM D638	1850 MPa
Elongation at Break	ASTM D638	5.60%
Shore Hardness	Scale D	86




OPTIMAL SCANNING

Models have smooth, hard surfaces for optimal scanning



ACCURATE & PRECISE

Highly accurate detail in printing



PAIRS WITH: PRO-FORM[®]

Check out our line of Pro-form thermoforming materials at dental.keystoneindustries.com

keystone
Industries



keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027

EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

REV 0 05/2025