



KeyOrtho IBT™

KeyOrtho IBT is ideal for printing orthodontic transfer trays for positioning brackets. KeyOrthoIBT's translucency allows for precise positioning before and during bracket placement. Its high tensile strength allows for easy application of brackets in a single step, reducing chair time to minutes. It is flexible for easy removal, and works with all general bracket systems. It is also odorless and tasteless for patient comfort.

#4200008 (0.5 kg bottle) | #4220008 (1 kg bottle)

PROPERTY	METHOD	KEYORTHO IBT
Tensile Strength	ASTM D638	8.0 - 10.5 MPa
Tensile Young's Modulus	ASTM D638	15.5 - 31.0 MPa
Elongation at Break	ASTM D638	>130%
Biocompatibility	ISO 10993-1	Pass
Viscosity @ 25°C	ASTM 7867	<1000 cP
Orthodontic Adhesive Release		Pass

FEATURES

- **BIOCOMPATIBLE MEDICAL DEVICE**
- **ODORLESS & TASTELESS**
- **ADHESIVES RELEASE EASILY FROM TRAY**
- **TRANSLUCENT FOR EASY POSITIONING**
- **PRECISE**
- **FLEXIBLE**
- **BPA, MMA & PHTHALATE FREE**

keystone[®]
industries



keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

REV 0 05/2025