



KeySplint Hard®

The newest member of the KeySplint family, KeySplint Hard is designed for rigid dental splints. KeySplint Hard provides therapy for bruxism and cases where tooth immobilization is needed.



BIOCOMPATIBLE



DURABLE



**EASY TO
POLISH**



**EASILY
CLEANED**



**ABRASION
RESISTANT**

FEATURES

BIOCOMPATIBLE MEDICAL DEVICE

- UNITED STATES FDA 510(K) CLEARED
- HEALTH CANADA LICENSED
- CE MARKED (CLASS IIA)
- HARD FOR RIGID SPLINTS
- BPA, MMA & PHTHALATE FREE

TINTED

#4200004 (0.5 kg bottle)
#4220004 (1 kg bottle)

CLEAR

#4200011 (0.5 kg bottle)
#4220011 (1 kg bottle)

PROPERTY	METHOD	KEYSPLINT HARD
Flexural Strength	ISO 20795-2	60-65 MPa
Flexural Modulus	ISO 20795-2	1510-1600 MPa
Water Sorption	ISO 20795-2	18 µg/mm ³
Elongation at Break	ASTM D638	9%
Biocompatibility	ISO 10993-1	Pass



keystone
industries

keyprint.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01

REV 0 05/2025